KP PLAST

Brilliant Red MG



◆Product summary

C.I.Vat Red 41
C.I.Solvent Red 242
C.I.Disperse Red 364
73300
Thioindigo dye
Powder
Bluish Red
100±2.5 (Relative ratio of absorbance for standard)
< 2.0 (Measurement by photo colorimeter with C/2 Light source)
< 1.0%
> 350℃
$0.10\sim0.30~\mathrm{g/cm}3$
< 2.0%
8.0 ~ 10.0


◆Solubility of Organic Solvents

Acetone	< 0.1
Methyl ethyl ketone	< 0.1
Ethanol	< 0.1
Ethyl acetate	< 0.1
Toluene	< 0.1 (g /100ml , 25℃±2℃)

◆Applicable resins

PS	0
AS	0
ABS	0
РММА	0
PC	0
PA	Δ
PET	Δ
PBT	Δ

◆ Visible light absorption spectrum



Notice: Please note that depending on your environment, there are some items and colors that can not be accurately displayed nor printed.